

TO REACH
HEALTHCARE PROFESSIONALS
WORLDWIDE...

ADVERTISE
ONLINE >

CLICK HERE

TO FIND
OUT MORE



IHE-ONLINE.COM

International Hospital Equipment & Solutions

The online source of medical device information connecting healthcare professionals to global vendors

IHE-online.com

[Magazine](#) | [This Week's News](#) | [E-products](#) | [Events](#) | [Contact us](#) | [My Ihe](#)



REKA Health System - Empowering Care



The new REKA Health System consists of ECG acquisition hardware, the E100 device, a PC or mobile OS app, and a cloud-based application platform ("REKA Health Cloud"). It is able to facilitate remote monitoring of patient's heart rhythms, by acquiring a 30-second record of the patient's ECG at the patient's discretion, after which the ECG record is uploaded to the REKA Health cloud using a PC or mobile application. These records can be instantly reviewed and diagnosed remotely by the patient's doctor, after which a response can be sent back to the patient via email or text message. REKA Health's technology offers several vertically integrated advantages for error elimination which is a key factor in remote patient monitoring. The E100 is a mobile pocket size ECG event recorder, designed to eliminate patient usage errors, and enables the reach to a larger user base. It works with laptops, PCs, and smart phones such as BlackBerry, Android and iPhone. For more information please visit the REKA Health website.

REKA Health
Haagwinde 20
5262 KZ Vught
The Netherlands
Tel.: +31-7365 76976
Fax: +31-7365 76977

